

# SWEET ACID FACE

## Regenerating cosmetic facelift

 Bottle 30 ml

Cod. 10103

### CHARACTERISTICS

Soluzione PH 3

### ACTIVE PRINCIPLES

**GLYCOLIC ACID**  
**FRUIT ACID**  
**KOJIC ACID**

### TREATMENT

#### **Smoothing, lifting, moisturizing, whitening, anti-wrinkle.**

Professional cosmetic facelift: the precious Alpha Hydroxy Acids help to reduce wrinkles, help to eliminate dead surface cells and rebalance the skin hydration.

**IMPORTANT:** in order to remove the ACIDO DOLCE VISO immediately after the end necessary, it is essential to prepare, before treatment, of gauze soaked in detergent solution.

Spread an even layer of ACIDO DOLCE VISO, being careful and avoiding that the product coli.

Start the application, starting with the less sensitive areas. Leave on for 30 seconds to 3 minutes, depending on skin sensitivity. Remove the acid with the detergent solution already prepared.

- IS RECOMMENDED USE OF LATEX GLOVES.
- FOR SENSITIVE, BEFORE, THE TEST PRODUCT ON A SMALL PORTION OF SKIN (AND ARM).
- IN CASE OF REACTION TO SUSPEND TREATMENT.
- TREATMENT EXCLUSIVE FOR BEAUTY CABINS.
- DO NOT APPLY ON EYE AND LIP AREA, ON INJURIES AND EXTERNAL MUCOSAE.

### INGREDIENTS

AQUA, GLYCOLIC ACID, GLYCERIN, PROPYLENE GLYCOL, SODIUM HYDROXYDE, KOJIC ACID, VACCINIUM MYRTILUS EXTR., CITRUS MEDICA LIMONUM EXTR., SOLANUM LYCOPERSICUM EXTR., CITRUS GRANDIS, MALIC ACID, CITRIC ACID, HYDROXYETHYLCELLULOSE, ETHYLPARABEN, METHYLPARABEN, PROPYLPARABEN, BUTHYLPARABEN, PHENOXYETHANOL, IMIDAZOLIDINYLU-REA.

